

LM 1

The in-lab mill for high productivity environments.



- Fast and precise 5-axis milling
- Broad range of indications with wet and dry milling capabilities
- Accepts various file formats and materials
- Ease of use with robotic 10-disc changer
- Open system integrated in the DTX Studio™ platform

KaVo LM 1. Technical specifications.

Wet/Dry milling	Wet and dry
Number of axes	5
Dimensions mm	576 width x 704 height x 581 depth
Machine weight	Approx. 140 kg
Blank size supported	98.5 mm, 10 to 40 mm height
Block size supported	40 mm max. length
Number of blocks supported	Max. 6
Automatic disc changer capacity	10 discs
Maximum spindle speed	80,000 rpm
Material support	Composites, LiSi2, Ti, PMMA, ZrO2, CoCr

Subject to technical modifications

Connecting to the DTX Studio™ platform.



For more information please register at: www.kavo.com

The products, equipment and services illustrated and described in this catalogue are not available in all countries.

All information conforms to knowledge at the time of printing. KaVo Dental GmbH accepts no liability for any deviation from the illustrations in terms of colour or shape, or any errors or print errors and reserves the right to make changes to the prospectus at any time. Full or partial reprinting is only permitted with permission from KaVo Dental GmbH.