


LS 3

The easy-to-use color desktop scanner you have been waiting for.



- High-speed color scanning, remarkable accuracy
- Touch screen display for ease of use
- Easy articulator scanning and SmartScan fixator recognition
- Connecting to the DTX Studio™ platform

KaVo LS 3. Technical specifications.

| | |
|-----------------------------------|------------------------------------|
| Dimensions* mm | 416 width x 412 height x 384 depth |
| Measurement field* (X x Y x Z) mm | 60 x 45 x 60 |
| Camera pixel | 2.8 MP |
| Accuracy* (ISO 12836) | Up to 4 µm |
| Measurement time complete jaw* | 19 sec. scanning |
| Measurement time single tooth* | 30 sec. scanning |
| Measurement time 3-unit bridge* | 43 sec. scanning |
| Weight | 35 kg |

Subject to technical modifications

* Ongoing testing against ISO standards pending

Connecting to the DTX Studio™ platform.



For more information please register at: www.kavo.com

The products, equipment and services illustrated and described in this catalogue are not available in all countries.

All information conforms to knowledge at the time of printing. KaVo Dental GmbH accepts no liability for any deviation from the illustrations in terms of colour or shape, or any errors or print errors and reserves the right to make changes to the prospectus at any time. Full or partial reprinting is only permitted with permission from KaVo Dental GmbH.